




KZ-700

Specifications

Type	KZ-700
Measurement method	Auscultatory measurement
Numerical display measurement range	Pressure: 0 - 300 mmHg Pulse: 30 - 200 beats / minute
Pressure bar display measurement range	Pressure: 20 - 280 mmHg
Measurement accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Power supply	2 x 1.5V alkaline batteries (LR6 or AA)
Upper arm circumference	23 - 33 cm using the medium cuff
Classification	Type BF 
EMC	IEC 60601-1-2: 2001
Operating condition	+10°C to +40°C / 30%RH to 85 %RH
Storage condition	-20°C to +60°C / 30%RH to 85 %RH
Dimensions	Approx. 98 (W) x 323 (H) x 66 (D)mm
Weight	Approx. 1.0kg - 101A

Options and accessories

AX-KZ-700-01	Medium cuff 23 cm to 33 cm
AX-KZ-700-02	Small cuff 16 cm to 23 cm
AX-KZ-700-03	Large cuff 33 cm to 45 cm
AX-KZ-700-10	Rubber bulb assy. (Rubber bulb, exhaust valve & end valve)

SUZUKEN COMPANY

8 Higashi Kataha-machi
Higashi-ku, Nagoya
Aichi 461-0015 JAPAN
TEL: +81-52-950-6327 FAX: +81-52-962-7440
<http://www.suzuken-kenz.com>