


Options & Accessories

◆ OS-30 Printer Unit

	2008 Y 9 M 18 D 15 H 10 m
Systolic	136 mmHg
Diastolic	83 mmHg
Map	99 mmHg
P.R.	74 b.p.m

Single result

Measurement Log



BP Trend Graph

Date	Time	Sys / Dia (Mean)	P.R.	Remark
10/28	17:17	126 / 77 (92)	63	
10/28	17:27	120 / 79 (91)	58	
10/28	17:28	125 / 83 (97)	58	
10/28	17:29	126 / 84 (98)	59	
10/28	17:30	119 / 80 (92)	58	
10/28	17:31	117 / 79 (90)	56	
10/28	17:32	118 / 75 (87)	58	
10/28	17:33	109 / 61 (76)	65	



The thermal printer unit is detachable from the back side of BPM OS-30.

Latex-free cuff Adult, Thigh type	38-50 cm
Latex-free cuff Adult, Large	31-40 cm
Latex-free cuff Adult, Medium	23-33 cm
Latex-free cuff Adult, Small	17-25 cm
Latex-free cuff Child, Medium	12-19 cm
Latex-free cuff Child, Small	8-13 cm
Extension Air Tube	250 cm



Bed-pipe Clamp HM-14



Bed-pipe Clamp HM-15



Mobile Cart HM-13

◆ Standard Accessories

Recording Paper MR-58 type 2 58 mm × 25 m	2 rolls
AC Power Cable	1 piece
Latex-free cuff Adult, Medium	1 piece
Extension Air Tube 100 cm	1 piece



Remote Switch

◆ Specifications

Power Supply	AC100-240 V 50/60Hz 30 VA	Surge Protection Class / Level	Class I B level equipment
Operation Mode	Single Measurement Interval Measurement	Compatible standard	Safety: JIS T 0601-1:2005 EMC: IEC060601-1-2:2004
BP Measurement Method	Oscillometric	Measurement Range	BP: 10-250 mmHg Pulse: 30-200 b.p.m
Display Range	BP: 0 or 10-299 mmHg	Sound	PCM composite
Transmission	Serial communication I/F x 1 Remote switch I/F x 1	Memory	100 measurements (displayed as lists or trend graphs)
Operational Temp.	Temp: 10-40°C Humidity: 30-85%	Temp. during Storage or Shipping	Temp: 10-60°C Humidity: 30-95%
Weight	Approx. 2.2 kg ±0.5 kg	Outer Dimensions	227(W) x 217(D) x 262(H) mm
Printer Unit Weight	Approx. 300 g	Printing Method	Line thermal

SUZUKEN COMPANY LIMITED

8 Higashikataha-machi
Higashi-ku, Nagoya
Aichi 461-0015 JAPAN
TEL: +81-52-950-6327 FAX: +81-52-962-7440
<http://www.suzuken-kenz.com>