

**Kenz BPM SP-1**  
**Automatic Blood Pressure Monitor**  
 More practical and easier operation  
 Superior accuracy of measurement

**Specifications**

<b>Measurement</b>	Oscillometric method
<b>Inflation pressure</b>	Automatic adjustment from 0 to 280 mmHg
<b>Display range</b>	Systolic 40 ~ 250 mmHg Diastolic 30 ~ 210 mmHg Pulse 30 ~ 180 mmHg
<b>Accuracy</b>	$\pm 4$ mmHg Pulse : $\pm 5$ %
<b>Deflation pressure</b>	Constant exhaust speed control
<b>Operation ambience</b>	+10 °C ~ +40°C
<b>Storage ambience</b>	-10 °C ~ +50°C
<b>Pressure detection</b>	Semiconductor pressure sensor
<b>Printer</b>	Thermal dot array
<b>Display</b>	LED digital display (Systolic, diastolic, mean, pulse, date/time)
<b>Power consumption</b>	32VA
<b>Power supply</b>	100~120V or 220-240V 50/60 Hz
<b>Dimensions</b>	450(W)x310(D)x280(H) mm
<b>Weight</b>	Approximately 9.0 kg
<b>Standard accessories</b>	Dust cover, Recording paper, AC power cable, operation manual